



OXYGEN PARTS INC.



Thank you for trusting Oxygen Parts Inc. with your oxygen needs!

We're here to help!

Feel free to contact us if you need help getting set up, are experiencing issues with your machine, or need help finding replacement parts or accessories.

ORDER PARTS



PHONE: 833-360-2155



WEBSITE: www.OxygenParts.com



EMAIL: CustomerService@OxygenParts.com



INTRODUCING Stratus 5

The 5 LPM Stationary Oxygen Concentrator from React Health

- Slim line design for a smaller footprint in the home
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet operation - compressor is specially placed, and padded to reduce noise
- Simple, intuitive operation



REACTHEALTH



Stratus 5L
Designed with the Patient in Mind

STRATUS 5L

Designed with your business in mind



Simple to service

Easy access to internal parts

3-year standard warranty

PRODUCT SPECIFICATIONS

Rated power	320W
Voltage	120V ± 312V, 60HZ ± 31HZ
Flow rate (L/min)	0.5 - 5.0
Oxygen Concentration	(93% ± 35%)
Sound pressure level dBA	(≤45)
Sound power level dBA	(≤55)
Net weight	34.1 lbs
Dimensions	15.3 x 9 x 23.6 inches



STRATUS 5 ALARMS

- Abnormal temperature alarm
- Power failure alarm
- Low flow alarm
- Low oxygen purity alarm
- Pressure alarm
- Thermistor fault alarm

PRODUCT CODES

Part Numbers

STR1005 - Stratus 5 Stationary Oxygen Concentrator

Replacement Part Numbers

STR1020 - Stratus 5 Replacement HEPA Paper Filter

STR1030 - Stratus 5 Replacement Sponge Filter

STR1040 - Stratus 5 Replacement Bacteria Filter

Oxygen Accessory Part Numbers

O2S2007 - 3B Standard 7' Oxygen Cannula

O2U2007 - 3B Ultra Soft 7' Oxygen Cannula

O2T2025 - 3B 25' Oxygen Tubing

O2C2030 - 3B Oxygen Tubing Connector



REACTHEALTH

info@reacthealth.com | www.reacthealth.com

React Health is formerly known as 3B Medical, Inc.