



OXYGEN PARTS INC.



Thank you for trusting Oxygen Parts Inc. with your oxygen needs!

We're here to help!

Feel free to contact us if you need help getting set up, are experiencing issues with your machine, or need help finding replacement parts or accessories.

ORDER PARTS



PHONE: 833-360-2155



WEBSITE: www.OxygenParts.com



EMAIL: CustomerService@OxygenParts.com



INTRODUCING
Stratus 5

All new 5lpm oxygen concentrator
from 3B™ Medical

- Slim line design for a less obtrusive profile
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet sieve bed exchange
- Designed for high reliability

STR1005 - Stratus 5 Oxygen Concentrator

 **3B™ Medical, Inc.**



Stratus 5



Easy Access

to internal components for servicing

Noise Suppression

Compressor is specially placed and padded to dampen noise levels

Simple

intuitive screen.



Model STR1005

Rated power W	320
Voltage	120V±12V,60HZ±1 HZ
Flow rate L/min	0.5~5
Concentration Rated flow	93%±5%
Sound pressure level dB A	≤45
Sound power level dB A	≤55
Net weight kg	15.5
Dimension in	15.3 × 9 × 23.6

FEATURES:

Abnormal temperature alarm; Low oxygen purity alarm; Power failure alarm; Timer; Pressure alarm; Low flow alarm; Thermistor fault alarm

When the nominal pressure of the oxygen output port is zero, the oxygen concentration is 93%±5%

Product Codes

Replacement Parts

STR1020 - Stratus 5 Replacement Hepa Paper Filter

STR1030 - Stratus 5 Replacement Sponge Filter

STR1040 - Stratus 5 Replacement Bacteria Filter



O₂ Accessories

O2S2007 - 3B Standard 7' Oxygen Cannula

O2U2007 - 3B Ultra Soft 7' Oxygen Cannula

O2T2025 - 3B 25' Oxygen Tubing

O2C2030 - 3B Oxygen Tubing Connector



3B Medical, Inc.

Info@3Bproducts.com www.3Bproducts.com