



OXYGEN PARTS INC.



Thank you for trusting Oxygen Parts Inc. with your oxygen needs!

We're here to help!

Feel free to contact us if you need help getting set up, are experiencing issues with your machine, or need help finding replacement parts or accessories.

ORDER PARTS



PHONE: 833-360-2155



WEBSITE: www.OxygenParts.com



EMAIL: CustomerService@OxygenParts.com

**Proven
best-in-class
reliability,
year after year.**



Oxlife INDEPENDENCE[®]

Portable Oxygen Concentrator

Our flagship device, the Oxlife **INDEPENDENCE**, sets the standard for **long operating life, reliability, durability and 24/7 performance.**

Equipped with DNA Technology, it is the **DME choice for hospital discharge** and the hospital-to-home transition, meeting the majority of long-term oxygen therapy requirements.



Continuous flow up to 3 LPM
and 12 pulse settings



Proven 24/7 use, year after year
backed by a five-year warranty




Use as POC & stationary
dual billing (E1390 & E1392)



No tanks to refill
or ongoing tank deliveries



24/7 remote monitoring
with DNA Technology




Long-life, dual battery
operation of up to six hours



Intuitive, easy to use
LCD display and setting controls





DNA TECHNOLOGY®

Remote asset & clinical monitoring

O₂ Concepts' proprietary, device-initiated connectivity solution

GPS-enabled inventory control
and loss prevention technology

24/7 remote monitoring
for device performance

Easily access data from any device
via secure, cloud-based platform

Proactively promote enhanced
patient therapy adherence

Reduce field calls with intuitive
customer support resources

Manage hospital DME closet
inventory with device Geofence



OXLIFE INDEPENDENCE® PRODUCT SPECIFICATIONS¹

Dimensions	Height: 20.2" • Width: 10.9" • Depth: 9.45"
Weight²	Device: 16.7 lbs
Oxygen Flow Rate (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)
Battery Duration³	5.75 hours @ 2 pulse at 20 bpm 2.5 hours @ 2 LPM continuous
Battery Recharge Time⁴	2.5 hours
Oxygen Concentration	87% - 95% at sea level at all settings
AC Power	Universal Input: 100 - 240V, 50/60Hz 50 - 60 Hz • 2.0A maximum
DC Power	12 - 15 VDC • 15 amps

¹ O₂ Concepts® reserves the right to revise specifications at any time without notice ² Production weight excludes cart kit
³ Battery duration presented for two-battery system ⁴ Battery charge time increases if device is operating while charging

Smarter, more efficient and more reliable oxygen therapy than ever before

A legacy of design innovations, cutting edge technology and unparalleled support

Relied upon by leading academic and research medical centers, and by hospital groups and homecare providers nationwide

The most cost-effective device & operating model

A simplified, effective hospital/clinic discharge process

The ability to remotely monitor and locate every device in your fleet

Designed for the reliability and durability needs of the DME industry

Contact us at **(877) 867-4008 (prompt 6)** or at **info@o2-concepts.com** for more information or to consult with one of our territory managers.

o2-concepts.com



O₂ Concepts®, Oxlife INDEPENDENCE®, DNA Technology® and Energy Smart Technology® are trademarks of O₂ Concepts®, LLC and may not be reproduced without express consent.
© 2024 O₂ Concepts®, LLC. All rights reserved.
900-1019 Rev A