



OXYGEN PARTS INC.



Thank you for trusting Oxygen Parts Inc. with your oxygen needs!

We're here to help!

Feel free to contact us if you need help getting set up, are experiencing issues with your machine, or need help finding replacement parts or accessories.

ORDER PARTS



PHONE: 833-360-2155



WEBSITE: www.OxygenParts.com



EMAIL: CustomerService@OxygenParts.com

Oxlife **INDEPENDENCE**[®]

PRODUCT SPECIFICATIONS¹

Dimensions	Height: 20.2" • Width: 10.9" • Depth: 9.45"
Weight²	Device: 16.7 lbs
Oxygen Flow Rate (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments
	Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)
Battery Duration³	5.75 hours @ 2 pulse at 20 bp
	2.5 hours @ 2 LPM continuous
Battery Recharge Time⁴	2.5 hours
Oxygen Concentration	87% - 95% at sea level at all settings
AC Power	Universal Input: 100 - 240V, 50/60Hz
	50 - 60 Hz • 2.0A maximum
DC Power	12 - 15 VDC • 15 amps
Operating Temperature	50°F (10°C) to 104°F (40°C)
Altitude Operating Range	0 - 13,123 feet (0 - 4,000 meters)
Transportation/Storage Humidity	0 - 95% non-condensing
Connectivity - DNA Technology	O ₂ Concepts [®] proprietary, device-initiated, asset and clinical management platform utilizing the Verizon [®] cellular network

¹ O₂ Concepts[®] reserves the right to revise specifications at anytime without notice

² Production weight excludes cart kit

³ Battery duration presented for two-battery system

⁴ Battery charge time increases if device is operating while charging



 **O₂ Concepts[®]**

O₂-Concepts.com | 877.867.4008

© 2022 O₂ Concepts, LLC; all rights reserved; O₂ Concepts[®] and Oxlife **INDEPENDENCE**[®] are registered trademarks of O₂ Concepts, LLC and may not be reproduced without express consent.