



OXYGEN PARTS INC.



Thank you for trusting Oxygen Parts Inc. with your oxygen needs!

We're here to help!

Feel free to contact us if you need help getting set up, are experiencing issues with your machine, or need help finding replacement parts or accessories.

ORDER PARTS



PHONE: 833-360-2155



WEBSITE: www.OxygenParts.com



EMAIL: CustomerService@OxygenParts.com

Oxy HOME[®]

Stationary Concentrator

Keep Going.sm

NEW

BENEFITS

- Quiet Operation
- Small Footprint
- Purity Check by Percentage
- Alarm Indicators: No Flow, Low/High Pressure, Purity, Power
- Twin Blower Cooling Design
- Contemporary Design
- 3-Year Replacement Warranty

PRODUCT SPECS **Part Number 1400-5100**

Flow Specifications	0 - 5 LPM
Oxygen Concentration	93% (+3%, -3%)
Oxygen Output Pressure	8 (+2%, -2%) psi
Rated Voltage/Frequency	120 VAC ±10% 60Hz
Power Consumption	350 Watts
Sound	≤ 45 dB(A)
Dimensions	13.4" W x 11.8" D x 25.6" H 34 cm x 30 cm x 65 cm
Weight	41.8 lbs (19 kg)



1400-5100

CONTACT US TODAY TO LEARN MORE

1-833-360-2155 • CustomerService@OxygenParts.com • www.OxygenParts.com