



## NIDEK NUVO 10

115 Volt  
60 Hz

230 Volt  
50 Hz

OCSI\* Part no.

595

1005

\* Oxygen Concentration Status Indicator

- ✓ Built on the reliable Nuvo platform
- ✓ Easy access to filter and hour meter
- ✓ Superior grade of molecular sieve
- ✓ High efficiency compressor

**Average power consumption** 700 Watts / 600 Watts

**Oxygen Concentration** 90% (+5% / -3%) at 10 LPM

**Liter Flow** 2 to 10 LPM

**Outlet Pressure** 20 PSIG (1.38 BAR)

**Weight** 54 lb / 24 kg

**Dimensions** 16"W x 16"D x 28"H  
39 cm x 40 cm x 71 cm