



OXYGEN PARTS INC.



Thank you for trusting Oxygen Parts Inc. with your oxygen needs!

We're here to help!

Feel free to contact us if you need help getting set up, are experiencing issues with your machine, or need help finding replacement parts or accessories.

ORDER PARTS



PHONE: 833-360-2155




WEBSITE: www.OxygenParts.com



EMAIL: CustomerService@OxygenParts.com

1025DS, 1025KS, 1025UK With Auxiliary Oxygen Port



PART DESCRIPTION	1025DS	1025KS	1025UK
1 Accumulator Tank	1025D-610	1025D-610	1025D-610
2 Auxiliary Oxygen Port	525DD-635	525DD-635	525DD-635
Cabinet Parts:			
3 Base	1025D-603	1025K-603	1025K-603
4 Compressor Box w/ Foam	1025D-615	1025D-615	1025D-615
5 Filter Door	525DD-639	525DD-639	525DD-639
6 Front Cover	525DD-611	525DD-611	525DD-611
7 Rear Cover	525DD-612	525DD-612	525DD-612
8 Cabinet Screws	525DD-636	525DD-636	525DD-636
9 Cable Tie (Large)	505DZ-617	505DZ-617	505DZ-617
10 Capacitor (Motor Start/Run)	1025DS-616	1025KS-616	1025UK-616
11 Caster, Non-locking	501DZ-603	501DZ-603	501DZ-603
12 Bed Check Valve	PVO2D-607	PVO2D-607	PVO2D-607
13 Circuit Breaker	1025D-613	1025K-613	1025K-613
14 Cooling Fan	1025D-634	1025D-634	1025D-634
15 Cooling Fan Guard	N/A	N/A	N/A
16 Compressor	1025D-625	1025K-625	1025K-625
17 Compressor Mounting Plate Package	1025D-632	1025D-632	1025D-632
18 Compressor Rebuild Kit	1025D-643	1025K-643	1025K-643
19 Exhaust Muffler	1025D-705	1025D-705	1025D-705
Filters:			
20 Cabinet Air Filter	303DZ-605	303DZ-605	303DZ-605
21 Intake Bacteria Filter	1025D-605	1025D-605	1025D-605
22 Compressor Filter	1025D-682	1025D-682	1025D-682
23 Final Bacteria Filter	PV5LD-651	PV5LD-651	PV5LD-651
Fittings:			
24 Accumulator Tank "Tee" (not shown)	444-582	444-582	444-582
25 Flow Meter	1025D-607	1025D-607	1025D-607
26 Purge Manifold (not shown)	N/A	N/A	N/A
27 Heat Exchanger Tube	1025D-655	1025D-655	1025D-655

PART DESCRIPTION	1025DS	1025KS	1025UK
Hoses:			
28 1/8" ID (Blue) (4' Lg)	444-554	444-554	444-554
29 Compressor Intake Hose (Silicone)	1025D-657	1025D-657	1025D-657
30 1/2" ID Silicone Braided (2' LG)	505DZ-634	505DZ-634	505DZ-634
Hose Clamps:			
31 Plastic (1/4" ID Hose) (not shown)	N/A	N/A	N/A
32 Ladder Clamp (1/2" ID Hoses)	444-566	444-566	444-566
33 IEC Connector	N/A	N/A	N/A
34 Motor Mounts	1025D-609	1025D-609	1025D-609
35 Oxygen Outlet Port	525DD-606	525DD-606	525DD-606
36 PC Board	1025D-622	1025D-622	1025D-622
37 Power Cord w/ Strain Relief	PV5LD-618	VARIOUS	VARIOUS
38 Power Cord Strap	MC29D-657	MC29D-657	MC29D-657
39 Power Cord Strain Relief	505DZ-645	N/A	N/A
40 Power Switch	1025D-508	1025D-508	1025D-508
41 Pressure Regulator (8.5 PSI)	MC29D-612	MC29D-612	MC29D-612
42 Pressure Regulator (20 PSI)	1025D-612	1025D-612	1025D-612
43 Pressure Relief Valve	1025D-614	1025D-614	1025D-614
44 Purge Harness Assembly	1025D-617	1025D-617	1025D-617
45 Rotary Valve	515ADZ-702	515ADZ-702	515ADZ-702
46 Sieve Bed	1025D-619	1025D-619	1025D-619
47 Top Cover Label	1025D-631	1025D-631	1025D-631
Wire Harness:			
48 Main Harness	1025D-623	1025K-623	1025K-623
49 Rotary Valve Harness	525D-621	525D-621	525D-621