

So much more than just
a medical device.

OxyHome® allows you
to enjoy life at home—
without interruption.

Oxy
HOME®

Oxy
HOME®

5L Stationary Concentrator

So quiet... you won't
even know it's on.

Contact Your Provider for More Information:

OXYGEN PARTS INC.
→ Take control

Oxygen Parts Inc.

6507 Jester Blvd

Suite 510-M

Austin, TX 78750

Hours: Mon-Fri 9AM-5PM CT

Toll-free: 833-360-2155

Austin: 512-372-1647

Fax: 800-418-6998

Email: help@oxygenparts.com

Website: www.oxygenparts.com

Connect With Us



oxygo.life

PN B150-0525



oxygo.life

oxygo.life

The Quiet Concentrator



The Future of Stationary Oxygen

The new OxyHome® Stationary Concentrator seamlessly blends style and functionality with its innovative design. Its compact size ensures that it fits effortlessly into any living space. But what truly sets it apart is its near-silent operation.

Enjoy daily activities without interruption; OxyHome is so quiet, you won't even notice it's running!

OxyHome Benefits

- Home Oxygen Filling Compatible
- Quiet operation
- Small footprint
- Alarm indicators: no flow, low pressure, purity, power
- Twin blower cooling design
- Contemporary look
- 3-year replacement warranty

LEARN MORE AT: oxygo.life



OxyHome Specifications

Part Number	1400-5100
Dimensions	13.4"(W) x 11.8"(D) x 25.6"(H) 34cm x 30cm x 65cm
Weight	41.8 lbs (19 kg)
Flow Spec.	0 - 5 LPM
Oxygen Concentration	93% (+3%, -3%)
Oxygen Output Pressure	8 (+2, -2) psi
Rated Voltage/ Frequency	120 VAC ± 10% 60Hz
Power Consumption	350 Watts
Sound	≤ 45 dB(A)

- ✓ Enjoy uninterrupted, peaceful sleep without disturbances from a noisy machine
- ✓ Move freely around your home, engaging in daily activities
- ✓ Interact with family and friends without the constant background noise from your oxygen concentrator