



Model No. 1400-2000NG

Every Unit Comes With:

- Single battery
- Nasal cannula
- AC charger for home
- DC charger for car
- Carry bag with padded adjustable strap

Optional Accessories:

- Intermediate Battery
Model No. 1400-2020NG-4
- Double Battery
Model No. 1400-2010-8
- Desktop Battery Charger
Model No. 1400-2030

Getting your OxyGo FIT+™ is easy.

Simply Contact Your Provider:



6507 Jester Blvd Suite 510-M
Austin, TX 78750

Website:
www.OxygenParts.com

Email:
help@OxygenParts.com

Toll-free: 1-833-360-2155
Austin, TX: 512-372-1647



More Flow. Same Freedom.

Breathe easy with 4 pulse dose settings — customized to keep you moving.

Connect With Us



www.OxygenParts.com

Lighter. Smarter.

Built to Keep You Going.

Keep Going.sm



Now With **4** Settings!

Introducing The All New

OxyGO FIT+™

Compact, quiet, and made for patients who need reliable oxygen without compromise.

Weighing just 2.9 lbs with the standard battery, OxyGo FIT+™ is one of the lightest POCs with 4 flow settings — and now offers enhanced features for even more freedom. Whether you're out running errands or enjoying life on the move, the FIT+ helps you **Keep Going**.

Like all OxyGo® portable oxygen concentrators, the FIT+ creates medical-grade oxygen from the air around you — no tanks needed. And it's designed for 24/7 use, at home or away.

- Easy-to-read digital display
- Simple, intuitive controls
- Designed for 24/7 use
- Meets FAA-acceptance criteria for use on board aircraft
- Optional Bluetooth connectivity (U.S. and Canada)



Model No. 1400-2000NG

LEARN MORE: www.OxygenParts.com

Size, Weight	With Single Battery: (L) 5.8" x (W) 2.7" x (H) 7.5" 2.9 lbs. / 1.3 kg. With Double Battery: (L) 5.8" x (W) 2.7" x (H) 8.1" 3.4 lbs. / 1.5 kg.
Oxygen Flow	Pulse Dose Settings: 1 (210ml/min) 2 (420ml/min) 3 (630ml/min) 4 (840ml/min)
Power	AC Power Supply: 100-240 V, 50-60 Hz (auto-sensing for use worldwide) DC Power: Car / Boat / RV compatible
Recharge Times	Single Battery: Up to 2.7 hours Double Battery: Up to 3 hours, 18 minutes
Battery Life	Single Battery: Up to 3 hours Intermediate Battery: Up to 4.25 hours* Double Battery: Up to 5.75 hours*
Noise	39 dBa (at setting 2) Quiet and discreet!
Warranty	3-year lifetime warranty See OxyGo Warranty Statement for details.

* At flow setting 1 ** When used under 5 hours/day