



Model No. 1400-3000NG

Every Unit Comes With:

- Single battery
- Nasal cannula
- AC charger for home
- DC charger for car
- Carry bag with padded adjustable strap

Optional Accessories:

- Double Battery:
Up to 12.75 hours of power
Model No. 1400-3010-16
- Desktop Battery Charger
Model No. 1400-3030
- Backpack
Model No. 1170-3420
- Custom Cart
Model No. 1170-3442

Getting your OxyGo NEXT® is easy.

Simply Contact Your Provider:



6507 Jester Blvd Suite
510-M
Austin, TX 78750

Website:
www.OxygenParts.com

Email:
help@OxygenParts.com

Toll-free: 1-833-360-2155
Austin, TX: 512-372-1647

Introducing the NEXT generation



The new OxyGo NEXT® keeps delivering — so you can **Keep Going.**sm

Connect With Us



www.OxygenParts.com

Keep Going.sm



12+
Hour
Battery!*

Breathe Easy. Live Fully.



The **NEW** OxyGo NEXT® is designed to keep up with your life – not slow it down.

With six intelligent pulse dose settings and extended battery options, the new OxyGo NEXT delivers powerful oxygen therapy in a compact, lightweight design. Quiet, reliable, and easy to use, it gives you the freedom to keep moving, stay connected, and **Keep Going**.

OxyGo® portable oxygen concentrators (POCs) are lightweight medical devices that make doctor-prescribed medical oxygen from the air around you. You never need to worry about running out of oxygen or carrying heavy tanks again.

- Simple controls, bright digital display
- User-replaceable batteries & components
- Designed for 24/7 use, at home or away
- Meets FAA-acceptance criteria for use onboard aircraft
- Optional Bluetooth connectivity (U.S. and Canada)

Model No. 1400-3000NG

LEARN MORE: www.OxygenParts.com

Size, Weight	<p>With Single Battery: (L) 7.18" x (W) 3.27" x (H) 8.14" 4.8 lbs. / 2.2 kg.</p> <p>With Double Battery: (L) 7.18" x (W) 3.27" x (H) 9.02" 5.8 lbs. / 2.6 kg.</p>
Oxygen Flow	<p>Pulse Dose Settings:</p> <ol style="list-style-type: none"> 1 (210ml/min) 2 (420ml/min) 3 (630ml/min) 4 (840ml/min) 5 (1,050ml/min) 6 (1,260ml/min)
Power	<p>AC Power Supply: 100-240 V, 50-60 Hz (auto-sensing for use worldwide)</p> <p>DC Power: Car / Boat / RV compatible</p>
Recharge Times	<p>Single Battery: Up to 3 hours</p> <p>Double Battery: Up to 4 hours</p>
Battery Life	<p>Single Battery: Up to 6.5 hours</p> <p>Double Battery: Up to 12.75 hours*</p>
Noise	<p>39 dBa (at setting 1) Quieter than a library!</p>
Warranty	<p>5-year lifetime warranty See OxyGo Warranty Statement for details.</p>

* Double battery at flow setting 1